

Dosimetric Result of Magnetic Resonance/Computed Tomography Compatible Hybrid **Phantom for Magnetic Resonance Guided Radiotherapy**

Min-Joo Kim¹, Seu-Ran Lee² and Tae Suk Suh²

¹⁾Department of Radiation Oncology, Yonsei Cancer Center, Yonsei University College of Medicine, Yonsei University Health System, Seoul, Korea²⁾Department of Biomedical Engineering, Research Institute of Biomedical Engineering, The Catholic University of Korea College of Medicine, Seoul, Korea

Hypothesis & Purpose

Hypothesis and Introduction

• MRI-only based radiotherapy which means that application of magnetic resonance image into radiation treatment target delineation, has been increased because MR image offers superior soft tissue contrast compared to CT image, especially, for head and neck tumors, prostate, lung tumors and brain lesions • Hypothesize that one single phantom named hybrid MR/CT compatible phantom used for acquisition of MR and conventional CT image could be applied for MR-only radiation treatment in terms of target delineation and radiation dose calculation • Time consuming, extra costs and registration accuracy associated with multiple imaging modality remains to be solved to fully establish MR image based radiation treatment and another more critical challenge is lack of ED information of MR image for dose calculation

*****Dosimetric evaluation

- CT image as reference data for dose calculation
- Converted MR image set as a testing image set by using commercialized treatment planning system and converting to CT-like image
- Collapsed cone algorithm (Version 3.2) in RayStation was performed by applying 500 cGy on isocenter of 6 MV photon for SSD (source to surface distance) 90 cm and 30 x 30 cm² fields with 0° of gantry angle

- Develop various tissue equivalent samples on each MR and CT image and inserted in the hybrid MR/CT compatible phantom
- Represent could represent tissue equivalently measured data for each signal intensity (SI) on MR image and HU on CT image from region of interest from one single sample on MR and CT image
- Directly applied into the conversion process from HU to relative electron density in TPS for dose calculation using the relation between SI and HU from the same sample image section

* Purpose

Develop the hybrid MR/CT compatible phantom using tissue-equivalent materials on each MR and CT image for magnetic resonance only-guided radiotherapy

• The center coordinate of phantom external region of interest as iso-center

RESULTS & DISCUSSION



Fig.1 Geometry of the developed phantom: Total height and external diameter was decided by internal size of 8-channel MRI head coil. Two disks at both end and one disk in the middle of phantom can provide two different interspaces. 7 plugs can be located at one of the intersp ace between the disks. Thus, total 14 plugs can be inserted into phant om. The developed phantom was customized to fix 14 plugs without water leakage



Fig.2 14 different developed tissue-equivalent materials (Upper) and its MR image for each plug listed in Table 1 with TR 1500 ms, TE 35 ms (Lower)

MATERIALS & METHODS

Requirement for "MR/CT compatible phantom"

- Generally, production of MR phantom should be considerate carefully compared to making of CT phantom because the principle of MR image acquisition is much more complex than those of CT image acquisition and most of MR phantom was filled with water or copper sulfate for image quality control or MR spectroscopy
- Requirement for "MR/CT compatible phantom" was established throughout paper research

Chemical component for tissue equivalent material

- Carrageenan as gelling agent which was blended from various seaweeds to produce rigid gel and more elastic and resistant to cracking than agar gel, and enable for the creation of large phantom for over extended time
- Agaros as T2 modifier
- GdCl3 as T1 modifier
- NaN3 as an antiseptic which is highly soluble in water and very acutely toxic NaCl and deionized and distilled water



Fig.3 Dose difference from CT and MR image based dose distribution. Two slice image was selected and dose scale on figure was selected from -5, -4, -3, -2, -1, 0, 1, 2, 3, 4, 5% for maximum dose, relatively.

- The importance of MR image phantom for radiation oncology field has been suggested. Also, to the best of our knowledge, another possible application of those tissue-equivalent materials in the field of radiation oncology, especially for radiation dose calculation, has not been reported.
- ✤ A hybrid MR/CT compatible phantom for image acquisition of MR and CT was designed and investigated the relation between MR and CT image for MR-only based radiotherapy in terms of target delineation and radiation dose calculation. In addition, various tissue-equivalent materials for both MR and CT image to be inserted into the developed phantom were described in this study.
- The modulation method of CT image and HU value was successfully developed by using K_2CO_3 since the formula to quantitatively calculate the amount of K_2CO_3 in tissue-equivalent materials for simulation of each organ and control of the HU value was deducted and the non-effectiveness of K₂CO₃ for MR image was demonstrated.
- However, T1- and T2- relaxation time of each tissue-equivalent material could not be achieved. The

K2Co3 as electron density modifier

• Various mixture of chemical component to simulate human tissue on both MR and

CT image was tested by measuring T1, T2 relaxation time and signal intensity on MR image and HU on CT image and each value was compared with reference data T1-T2 measurement data was compared with reference data and the average of percentage difference was 2.86% and -13.37% for T1- and T2 relaxation time, respectively. With the 5% of tolerance for percentage difference, 10 measurement data among 24 was acceptable and tissueequivalent materials of liver was acceptable as tissue-equivalent materials.

2020 KSMP Fall Poster 0055 smiling.mjkim@gmail.com / +82) 10-26746091